Jateador  
Shot Blaster  
Jateador

# professional jet

**EDG**  
EQUIPAMENTOS

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Serial Number

## 1 Receipting


Dear customer, first of all we would like to thank you for acquiring EDG equipment, which are developed with the highest technologies and quality standards. However, to keep these characteristics, we would also like to warn you that some care should be taken when receipting the equipment. After receipting the equipment, check the general conditions of the package, if any damage is noticed contact the company immediately which will provide guidance on how to proceed. We remind you that the equipment is shipped at the purchaser's account and risk and it is insured by the conveyor.

Your Professional Jet package should include:


- a) The equipment itself.
- b) An activation pedal.
- c) Instruction manual.

## 2 Description

- 1) Abrasive materials selector;
  - 2) Air entry valve;
  - 3) 2mm movable air blaster;
  - 4) 1mm movable air blaster;
  - 5) x;
  - 6) Light bulb shield;

- 7) Manometer;
  - 8) Acrylic shield;
  - 9) X;
  - 10) Gloves;


## Notes

- Pressures over 5kg/cm<sup>2</sup> cause excessive and unnecessary wear to the equipment components.
- Clean the aspirator filter from time to time according to the usage frequency.
- To remove the material excess already used, remove the cap located at the bottom of the equipment. The inner grids can be easily removed for cleaning using a spatula.
- The light bulb used (pic.1-7) is the 9 watts without reactor.
- The movable air blast nozzles (pics.5-22, 23) are interchangeable and can be detached from their supports (pic.5-21). When reconnecting the nozzles, clean the screw avoiding the abrasive material to jam it. The nozzles inner diameters are: 1mm e 2mm.
- If a nozzle does not work:
- Tap it with a finger and activate the pedal.
- The used material has too big granulometry to the nozzle. Change the 1mm nozzle toe the 2mm one.
- The material is wet and there is water accumulation in the compressor and its piping.


## SPECIFICATIONS

Nominal supplies sources: 110 or 220Volts according to the ordering.  
 Consumption: 50 watts 50/60 Hz. (except the aspirator)  
 Maximum copressed air pressure: 06 kg/cm<sup>2</sup>, 90 lb/in<sup>2</sup>  
 Working pressure: 01 to 04 kg/cm<sup>2</sup>, 15 to 60 lb/in<sup>2</sup>  
 Maximum outflow: 100 litrers/minute


- 11) Tank "2";  
 12) Tank "1";  
 13) Air outflow for the aspirator (suction hose connection);

- 14) Power cable connection;  
 15) Activation pedal connector receptacle;  
 16) Main switch;  
 17) Aspiration speed activation button;  
 18)Aspirator electric socket;  
 19) Compressed air inlet;

- 20) Movable nozzles supporting cables;  
 21) 2mm movable nozzle;  
 22) 1mm movable nozzle;

### 3 Installing

- Professional Jet needs 4 kg/cm<sup>2</sup> or 60lb/in<sup>2</sup> pressure compressed air with 100 liters per minute yield.

Usually, a compressor with ¾ HP motor is enough.

- It is advisable to install a filter for eliminating water from the compressor outlet. Adjust the compressor to turn on with 60lb/in<sup>2</sup> and off with com 90lb/in<sup>2</sup>.

- Connect the compressor to the equipment using a proper PT 250 ¼ inches braided hose. (pic.4-20).

- WE remind you that for a good performance and little charge loss\*, the distance between the compressor and the air blast should not be longer than 5 meters. For longer distances the hose diameter should be widened to 3/8".

\* Air leakages can be understood as unnecessary energy waste, thus perform all connections in a way to avoid leakages.

- Load the left side (pic.3-13) "tank 1" with smaller granulometry materials and load the right side (pic.3-12) "tank 2" with bigger granulometry materials, do not let the material level to surpass the tank inner feeding tube.

- Next, close the compartment until feeling a slight resistance, but never tighten too much.

\* The system is self-sealing

- Avoid that abrasive materials get in contact with the screw and the lids.

\* If you wish, you can put some aluminum oxide or another mixture already used in the chamber to test the fixed nozzle (pic.1-6).

- Connect the activation pedal to its receptacle (pic.5-16).

- Connect the aspirator socket to the female plug located at the equipment backside (pic.4-19);

- Connect the aspirator hose to the air blast using the connection (pic.3-14).

- Connect the air blast plug to an electrical network compatible with the tension shown on the label (pic.4-15). It is advisable to follow the safety grounding norms for all electric equipment, to do so use a grounded rod or you building grounding. Never use the network neutral as grounding.

**Notes:** The aspirator is not provided with the equipment. Any vacuum cleaner should be appropriate to the process. The aspirator should have the same tension as the air blast working tension.

\* Charge loss is the difference between the working air blast pressure and the static air blast pressure. It happens because of the air passage resistance through the hoses, connections, etc.

### 4 Operation

- 1) Turn the aspiration speed adjustment button (pic.4-18) to its minimum;
- 2) Turn on the main switch (pic.4-17);
- 3) Turn the aspirator switch on;
- 4) Close the acrylic shield (pic.2-9);
- 5) Adjust the aspiration speed in a way that the gloves (pic.2-11) are slightly inflated on horizontal position.
- 6) Put your hands in the gloves and adjust the compressed air pressure according to the work to be performed turning the pressure switch (pic.1-1) and reading the value on the manometer (pic.2-8).
- 7) Select the abrasive material using the selector key (pic.1-2).
- 8) Activate the pedal directing the abrasive blast to the piece to be worked on.
- 9) To change materials, just turn the selector key (pic.1-2) to the desired position, "1" left tank, "2" right tank and "0" fixed nozzle for already used material.

**Notes:** Never direct the abrasive blast directly to the shield, because it will cause irreversible damages to the shield.

### 5 Maintenance

- Replacing the shield protection.

The Professional Jet shield is protected by a "Pac Film" plastic film, which can be easily found in supermarkets. After some time, the film gets opaque and it is necessary to replace it. Remove the old film, clean the shield with a wet cloth, put the new film on it and, using a stationery knife over the groove (pic.2-10), cut the exceeding material out.

- Keep the air entry valve (pic.1-3) with its latch inside the blast chamber.
- The gloves (pic.2-11) when worn out or torn can be easily replaced by stretching the wrist and adapting them to their support. The gloves used with the equipment are the big size yellow type.