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
--

## 1 Receiving

The equipment is packed separately from the vacuum pump, check if you have two boxes if you have acquired the furnace and the pump. Check the general conditions of the package, if any damage is noticed contact the shipping company. We remind you that the equipment is shipped at the purchaser's account and risk and it is insured by the shipping company.

**The FV Plus Automatic package should include:**

- a ) A furnace
- b ) A platform for ceramic burn;
- c ) A pack with 6 isothermal pins;
- d ) Instruction manual.

The vacuum pump package should include the vacuum pump and its accessories.

## 2 Main characteristics

- Quartz muffle with thermal insulation ceramic mat molded through low thermal mass vacuum, which insures perfect temperature homogeneity.
- 50 programs, being 49 available for customizing.
- Management by microprocessor.
- Easy operation and visualization of the burn parameters.
- Ambient working temperature at 1,100°C with maximum temperature limiter at 1,200°C
- Linear heating speed from 1 to 70.0°C/min.
- Thyristored power control
- Ceramic burn time from 0 to 30 minutes.
- Vacuum time from 0 to 30 minutes.
- Decreasing time visualization.
- Automatic or manual (which can be activated ant any stage of the process) vacuum operation
- Swinging tray, which does not let the pieces to be burnt vibrate.
- Programmable and independent closing and opening times.
- Safety systems protect the equipment from operation mistakes or defects.
- Flashing and beeping signs of all process steps.
- To avoid ambience heating and save energy, an alarm will sound after 2 minutes if the equipment is not being operated the muffle will be partially closed automatically.
- Muffle forced ventilation for quick cooling.
- Operation time totalizer.
- Low cost muffle replacing system.
- Electronic and heating systems forced ventilation keeps the components safe.

## 3 Installing

The numbers in parentheses refer to the picture that follows.

Your FV Plus Automatic should stay far from curtains and inflammable material. A furnace is a heat generator which needs to be dissipated; otherwise its components will be overheated. Therefore, place the equipment at a ventilated area which allows free air circulation. It is advisable to leave a 15cm clearance between the furnace and anything that can impair the ventilation. Place the furnace far from taps or basins which can sprinkle water over the equipment.

**Error 2** - Muffle elevator flaw during ascending procedure indicates mechanical problems. Check if there are no obstacles blocking the elevator movements.

**Error 3** - Elevator flaw during descending procedure indicates mechanical problems. Check if there are no obstacles blocking the elevator movements.

**Error 4** - heating speed equals zero. To fix the problem, type PROG and choose a value different from zero.

**Error 5** - Saving recipes out of the standard recipes for sintering/infiltration can not be saved on the positions reserved for ceramic recipes and vice-versa.

**Error 6** - Attempt to change the standard recipe 01. The standard recipe cannot be changed. Save the new recipe in another available position in the memory.

**Error 7** - The temperature has exceeded the maximum limit allowed when working with the vacuum (1,100°C) for this application. Check the programming, if the error persists call the technical assistance.

### 7.4 Important notes

FV Plus Automatic is calibrated to perform instant measurements of temperature..

-While the muffle is open or under closing procedure some temperature variations can be observed which are mainly caused by the ambient air in the furnace; this phenomenon is normal and does not interfere in the burn result.

- For burn processes without burn time, high heating velocities and large pieces, some temperature adjustments may be necessary for certain ceramic types.

- For a better result, it is convenient to use low thermal platforms and inert alumina pins of low thermal conductibility (isothermal pins).

#### Specification:

- Power supply: 110 or 220 volts. 50/60 Hz. Under request
- Maximum consumption: 1400 Watts.
- Fuses: 20 amperes to 110 volts, 15 amperes to 220 volts.
- Maximum operation temperature: 1100°C with the vacuum.

#### Equipment dimensions:

- |                            |                          |
|----------------------------|--------------------------|
| - Width.....               | 285mm                    |
| - Height when closed ..... | 320mm                    |
| - Height when opened ..... | 390mm                    |
| - Depth.....               | 405mm                    |
|                            | - Net weight..... 10.8Kg |

#### Package dimensions:

- |                         |                          |
|-------------------------|--------------------------|
| - 350mm x 420mm x 450mm | - Gross weight..... 12Kg |
|-------------------------|--------------------------|

For the burn and vacuum times course, the times can be visualized activating the Tempo de queima (Burn time) key (11) or Tempo de vácuo (vacuum time) key (3). The times will be visualized on the display decreasingly.

#### Resting position

After 2 minutes without any operations, the equipment gets into a resting position. The muffle is partially closed, and be kept that way until the next burn process. The display will show the "Prot" message. This process helps saving energy and avoids the unnecessary ambient heating.

To return to initial position use the Elevador manual (manual elevator) key up. If the operator does not the resting position procedure to happen, just put the muffle down manually a few centimeters using the Elevador Manual para baixo (manual elevator down) key (12)

## 7.2 Usage time

The furnace is equipped with a device that sums up the burn hours. From the moment the S/S key is activated, the hours used information is saved in an inner clock. This device is important to control the muffle useful life as well as the replacement parts. It also helps controlling energy waste etc.


To access the total usage time activate, on the parameters adjustment keyboard, the keys for thousands, hundreds and decimals following this sequence and then the Tempo de Vácuo

(Vacuum time) key, the value shown on the display is the hours total which is saved in the equipment memory. To exit this function, activate the Salvar (save) key. Before being shipped, your equipment is exhaustively tested; thus, values ranging from 1 to 10 working hours are common in new products.

## 7.3 Error messages


To avoid damage to the equipment and programming mistakes, several safety devices work in the flaws prevention system. Messages and error codes will be visualized on the display.

**Error 1-** Thermocouple flaw (temperature sensor). The display will show a code for the error and the process being performed is aborted. Contact the technical assistance.

Pic-1

- 1-1) Check if your network voltage is the same as the one shown on the furnace (2).
- 1-2) Install the furnace into an exclusive electrical network using 6mm wiring to a 110 volts network, or 4mm wiring to a 220 volts network.
- 1-3) Never connect the furnace to the same network in which other furnaces, compressors, electrical taps or any other high electricity consumption devices are already connected.

Pic-2

- 1- Fuse holder.
- 2- Network tension identifier/ power cable.
- 3- Communication for check-up (technical assistance)
- 4- Vacuum hose connection.
- 5- Vacuum pump socket.

**1-4)** Check if the outlet in which the furnace will be plugged in is in perfect condition and of good quality to bear at least 20 amperes.

**1-5)** Connect the grounding end (round pin) to a grounding rod and never to the electrical network neutral.

**Warning:** Not checking the procedures above will interfere on the equipment correct working, on its guarantee and also on the operator safety.

**1-6)** If your electrical network is very instable, which can be checked if the lights flash or there is intensity alteration, it is advisable to use a voltage regulator such as a saturated autotransformer with at least 2kw capacity. Never use, by any means, regulators used in computers.

**1-7)** Connect the vacuum hose to the furnace (4) and to the pump. Plug the vacuum pump power cable to the outlet which is at the backside of the furnace (5).

## 4 Operation

To better illustrate the operation, we bring a programming example for metaloceramic burn.

### Example:

Starting from the ambient temperature, go to T1 temperature at 450°C; furnace work entry temperature.

- Keep this temperature until the operator command to start the burn process.
- The muffle closing process should take about 5 minutes.
- The burn temperature is 935°C. (T2).
- The heating speed between T1 and T2 should be 55°C/minute.
- The burn time should be 5,35 minutes.
- The vacuum should be turned on during the muffle closing and be turned off 1.37 minutes after the burn time is started.
- After the process is over, the muffle should open in 1 minute.

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Choose the recipe number to be read using the parameters adjustment keys.

Activate the LER (READ) key (4). The recipe is ready to be used. To read a recipe it is necessary that it has been previously saved.

## 6.10 Burn

Pic-4

- So far, all parameters are programmed and saved:

- Wait until the temperature reaches the one programmed for T1. Aguarde a temperatura atingir aquela programada para T1.

- Place the pieces on the center of the burn platform (B), symmetrically leaving equidistant spaces between them for better heat distribution,

- Activate the start/stop key (14) once; the program will run.

- The muffle will be closed and the vacuum pump is activated. The vacuum is shown on the vacuometer and the temperature will rise until it reaches T2.

- After T2 is reached, the burn and vacuum times start counting; when it is done, the vacuum is released and the burn process continues without vacuum. When the burn time is over, the message “Ar” is shown and the muffle opens to its upper point.

- The equipment starts a cooling down process until the temperature reaches T1 again, so that it is ready for a new burn process.

- As the burn process is being performed, leds will indicate it.

- All parameters are saved in the memory. Similar and successive burn processes can be performed without the need of new programming.

- The equipment starts a cooling down process until the temperature reaches T1 again, so that it is ready for a new burn process.

- As the burn process is being performed, leds will indicate it.

## 7 Other functions

### 7.1 Burn and vacuum times visualization

## 6.7 Adjusting the muffle opening and closing times

The muffle closing and opening processes can be programmed independently. Activating the Abre/ Fecha (open/close) key (13) once the display will show the message "desc-1".

Activating the parameters adjustment key below "1" we choose the desired position according to the table below. In this case, on the position "5" the piece will take about 5 minutes to complete its entry.

Activating the Abre/ Fecha (open/close) key (13) twice the display will show the message "sob.d".

Activating the parameters adjustment key below "d" we choose the desired position according to the table

In this case, on the position "1" the piece entry will take about 1 minute to complete the opening course

### Muffle opening and closing modes table.

"0" = manual activating using the keys (12).

"D" = the muffle opens and closes non-stop.

On the next positions, the course times are approximately:

"1" = 1 minute

"2" = 2 minutes

"3" = 3 minutes

"4" = 4 minutes

"5" = 5 minutes

## 6.8 Saving the recipe

When the programming process is over we will save the recipe. As the recipe above is a sintering/infiltration one, we should choose the recipe number from 02 to 50.

Activating the PROG. Key (4), the display will show the message "C-01" which is the number of the standard recipe that comes from the factory.


Using the parameters adjustment keys we give a number to the recipe, for example, 02 and activate the SALVAR (save) key (4).

The recipe will be saved in the position "C-02" is the memory.

All parameters can be altered during any stage of the process except the burn and vacuum times when they are being performed.

		Muffle				
T1°C	T2°C	VEL.°C/m in	Burn Time min.	Vacuum time min.	Descent (Close)	Ascent (Open)
450	935	55	5,35	1,37	5	1

## 5 Programming

- (1) Parameters adjustment keys
- (2) Display.
- (3) Vacuum time key.
- (4) Program system entry key, save and read recipes.
- (5) Functions indicative leds.
- (6) Normal vacuum key.
- (7) Vacuometer
- (8) Temperature key 1.
- (9) Temperature key 2.
- (10) Heating speed key.
- (11) Burn time key (baseline).
- (12) Muffle manual opening and closing keys.
- (13) Muffle opening and closing programming keys.
- (14) Start / stop, key to begin and finish the burn process.
- (15) Equipment main switch.

## 6 Step-by-step description

### 6.1 Turning the FV Plus Automatic on

Turn the main switch located on the panel on (15); the display (2) will show the program used; next the message "Ar" will be shown for about 30 seconds while the equipment performs its checking-up; After that, the ambient temperature will be shown; if the muffle is close it will automatically open. The heating will start until it reaches the T1 temperature of the standard program that comes from the factory (400°C). The basic standard program which comes from the factory is saved in the position 01 of the memory; all programming will be done by altering its parameters and saving the new ones in other positions of the memory.

**Note:** Four seconds after the functions keys are activated and the parameters values are adjusted, the system will automatically return for the insertion of a new burn parameter.

When activating the PROG. (4) key, the display (2) will show the message "C-01", which is the number of the sintering/infiltration standard recipe that serves as a base for programming the next recipes, and then press the "LER" (4) (read) key.

### 6.2 T1 and T2 Temperature adjustment

- Activate the T1 (8), key and the display will show the programmed temperature of the standard recipe that comes from the factory "0400"
- Next, activate the parameters adjustment keys (4), one for each digit showing on the display the desired value of 0450°C. The furnace will start heating until it reaches this temperature..
- Activate the T2 (9) key and the display will show the programmed temperature of the standard recipe that comes from the factory "0960".
- Next, activate the parameters adjustment keys (4) one for each digit and set the desired value for T2 of 0935°C.

### 6.3 Adjusting the Heating Speed

- Activating the VELOC. (10) key the message "A060" °c/minute will show on the display. Which is the heating velocity for the standard recipe that comes from the factory.

- Activate, on the parameters adjustment keys (4), one for each digit and set the desired value of "A055" °c/minute.

### 6.4 Adjusting burn time

**Note:** The burn times for ceramic are indicated in minutes and seconds. Ex: 00.30 corresponds to thirty seconds; 05.35 correspond to five minutes thirty-five seconds.

- Activate the burn time key (11) once, the display will show "01.00" minute, which is the burn time for the standard recipe that comes from the factory
- Activate, on the parameters adjustment keys (4), one for each digit and set the desired value of "05.35" minutes

### 6.5 Vacuum time adjustment

- Activate the Tempo de vácuo (3) (vacuum time) key, and the display will show "00.00" minutes which is the standard recipe vacuum time that comes from the factory.
- Next, activate parameters adjustment keys (4), one for each digit showing on the display the desired value of "01.37" minutes.

**Note:** - The vacuum time can never be higher than the burn time, if the value for the vacuum time is higher than the value for the burn time, the furnace will determine for the vacuum time the same value as the burn time.  
- Vacuum time the same as "00.00" (zero) indicates that the vacuum will be performed on the heating ramp between T1 and T2. When the T2 temperature is reached the vacuum will be turned off.

### 6.6 Manual vacuum adjustment sim/não (yes/no)

The vácuo manual (manual vacuum) key (6) has two functions:

- Activate the manual vacuum key (6) when you want to interrupt the vacuum at any stage of the process, the pump will be turned off and the vacuum is undone.
- When a burn without vacuum is desired (oxidation) the vacuum time should be defined as "00.00" and the vácuo manual (manual vacuum) key (6) should be activated in a way that the led above it (5) is off.