*Solda plasma pontual micropulsada*  
*Micropulse weld plasma*  
*Solda de plasma pontual micropulsada*

 **MicroMelt**

**EDG**  
EQUIPAMENTOS

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|  |
|--|
|  |
|--|

## Presentation

In a lab, the blowtorch welding is a critical moment in a prostheses manufacturing process, it is inaccurate and changes the metal physicochemical characteristics creating weak points in structures, especially in implantology in which the welding is largely used. Besides, the blowtorches make use of gases which are dangerous for indoor usage.

Micromelt uses the welding process through micropulse plasma in which an alternated current creates plasma between a tungsten electrode and the piece surrounded by argon atmosphere. The plasma is accurate and reaches 3,000°C. By controlling the plasma intensity and duration we have a perfect heat source for accurate welds using less heating on the areas surrounding the welding pinch, keeping the base metal original characteristics. It can be used to bond all kinds of metal, including titanium and its alloys. As it does not use high tension and inflammable gases it is totally safe, and its adjustable parameters allow the repetition of the results no matter the operator.

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## 1. Reception

When receipting the equipment, before signing the reception counterfoil check the general conditions of the package, if any damage is noticed contact the company immediately which will provide guidance on how to proceed. We remind you that the equipment is shipped at the purchaser's account and risk and it is insured by the conveyor.

After opening the package, check its content:

- Equipment.
  - 2 grounding wires with clips
  - 4 tungsten electrodes
  - 1 hose for argon inlet.
  - 1 compressed air hose
  - 4 clamps.
  - 1 pedal with connection cable
  - 1 ceramic socket
  - 2 hand rests
- Instruction manual.

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**Warning:** Reading this manual carefully before operating the equipment is an owner responsibility. Any utilization different from what the equipment was designed for will be considered inappropriate usage, as well as not following the operation, safety and maintenance rules indicated by this manual.

## 2. Safety Recommendations

- 2.1.** Do not look at the plasma without the protection devices supplied with the equipment.
- 2.2.** Instruct all operators and assistants about the safety procedures.
- 2.3.** Never let lay or curious people operate the equipment, especially using it differently from what it was designed for.
- 2.4.** It is recommended the usage of protection gloves against possible cast metal or slag sparkling.
- 2.5.** Remember that you are operating electric powered equipment, so do not operate the equipment when barefoot or on wet floor because spurious currents can be formed. Prefer wearing rubber shoes.
- 2.6.** Check the grounding in which the equipment is connected. Make use of a professional electrician. The best ground resistance is 5 Ohms.
- 2.7.** The argon cylinder contains gas under high pressure, so follow your gas supplier safety recommendations.
- 2.8.** Do not let inside or next to the equipment containers, cloths or pieces of cotton which have inflammable liquid because weld spatters can set fire on them.
- 2.9.** Never weld containers which have or already had fuel in them, there are explosion risks.
- 2.10.** Pacemaker bearers should contact the pacemaker manufacturer for possible limitations.

## 8 Specification

- Nominal supply source: 220 V
- Quiescent current: 1,5 Amp.
- Pulse current: Max 25 Amp.
- Shielding gas: Argon.
- Minimum purity rate recommended: 99.99%
- Best purity rate: 99.999%
- Gas consumption: 5 liters/min.
- Electrode: Tungsten.
- Maximum compressed air pressure: 2atm ( 30 psig.)
- Dimensions: Width:0,42m ( 16,5")  
Height:0,29m ( 11,5")  
Depth:0,48m ( 19,0")  
Net weight:19 kg (43lb)  
Gross weight: 23 kg (51lb)

## 4.3 Argon installation

Micromelt needs argon for working, which is inert, non-flammable and protects the weld area against oxidation of the metals to be welded. Pic:3

Argon is supplied in pressurized cylinders with different capacities. A 3m<sup>3</sup> (100ft<sup>3</sup>) cylinder is the most indicated and enough for hundreds of welds.

Argon is supplied in different purity rates. Use preferably the 4.5 rate (99.995 % purity) or rate 5.0 (99.999 purity).

The argon cylinder should be attached to Micromelt through a specific pressure and outflow regulator to use with argon (Pic. 3) which can be obtained with your gas supplier.

The right gauge shows the gas pressure inside the cylinder, the left gauge shows the gas outflow which is regulated through the adjust knob.

Install the regulator to the cylinder and the regulator outlet to the argon equipment inlet pic. 2, to do so use the hose and clamps

Pic. 3

## 4.4 Compressed air

The Micromelt working does not depend on compressed air installation.

Compressed air is used for refrigeration of the pieces after welding when needed.

To install it connect the air hose to the inlet indicated in pic. 2. The pressure should be adjusted to 2Kg/cm<sup>2</sup> (30 psi) at most, for that use an air pressure regulator.

## 4.5 Pedal

The welding activation is done through a pedal connected to the back connector Pic. 2. The pedal has two functions: When it is pushed it discharges the argon flow and when the pedal is released it activated the plasma.

## 4.6 Tilting cabinet

Pic. 4

To facilitate the access to inner areas and for maintenance, the Micromelt cabinet opens up.


## 4.7 Torch

Micromelt is shipped with the torch assembled and ready to use. Pic. 5 shows the assembling sequence.

Tungsten electrode (1) ; Ceramic socket (2) ; Electrode bracket (3) ; Torch body (4) ; Mandrel (5) ; Back body (6)

- Screw the electrode bracket (3) to the torch body (4).
- Screw the ceramic socket (2).
- Insert the mandrel (5).
- Screw the back body without tightening (6).
- Insert the electrode (1) through the ceramic socket leaving a 10mm to 15mm (3/8" to 5/8") clearance beyond the socket.
- Tighten the back body firmly until the electrode is fixed.

## 6 Maintenance

Your Micromelt maintenance is very simple. You just need to keep the tungsten electrodes sharp for obtaining optimal welding quality.

Check the argon amount in the cylinder from time to time.

Keep the working area always clean and, for obvious reasons, never use inflammable products for cleaning.

After the usage, turn the equipment off and close the argon valve, because even small leakages after a long time can empty the cylinder charge.

Clean the magnifying lens and the photo sensitive crystal from time to time with a dry cloth.

## 7 Problems and solutions

| Problem   | Cause  | Solution  |
|---|--|---|
| - the equipment, even being with the main switch on, does not work. | -Poor contact in the electrical outlet.<br>- Fuse.                                 | - Check if the outlet is well connected and if the electrical network is energized.<br>- Check the fuse |
| - There is not plasma formation.                                    | - The grounding clip is not connected or the wire is split.                        | - Connect the grounding.<br>- Replace the wire.   |
| - Small and insufficient weld pinch.                                | - Low pulse value.   | - Turn the pulse value up.  |
| - There is plasma formation but no welding.                         | - Pointless electrode.<br>- The argon flow is not open                             | - Sharpen the electrode properly.<br>- Release the argon  |
| -There are fumes and too much oxidation.                            | - The argon flow is not open or its outflow is too low.                            | - Adjust the argon outflow.   |
| - There is dark oxide formation                                     | - Argon outflow is too high.   | - Adjust the argon outflow  |
| - The electrode sticks too easily.                                  | - Pointless electrode.<br>- There is some pressure from the piece on the electrode | - Sharpen the electrode.<br>- Touch the piece with the electrode without pushing                        |
| - The weld pinches are irregular.                                   | - The electrode is loose   | - Tighten the back body properly (Pic. 5).  |
| - The piece gets melted but there is not welding                    | - Pulse and depth adjustment too high for the piece size                           | - Adjust the pulse and depth values properly profundidade.  |
| - The pulse and depth values cannot be altered.                     | - The pedal is pushed, the values can only be altered when the pedal is resting.   | - Release the pedal when making any adjustments.  |

and the argon will keep flowing for some seconds.

- It is possible that on the first attempts the electrode gets bonded to the piece making it necessary to resharpen the electrode, make several attempts until you get practical,, remember that the secret is to touch the piece on the electrode without pressing it.
- Then put two pieces that have not been used together.
- Vary the pulse and the depth to feel what values are more appropriate.
- As Micromelt welds several different metals in the most varied situations, we will nt specify the pulse and depth parameter adjustments. With the practice, the operator will easily find the most appropriate values for each situation.


**Warning:** Micromelt is prepared for welding some large pieces as implant protocols, to do so it can be necessary using some adjustments close to the equipment maximum capacity (pulse 60, depth 20), under these conditions high energy plasma will be created generating a great amount of located heat. In this case, keep your fingers far from the welding pinch.

### 5.3 Some suggestions

- Clean the area to be welded.
- Always fix the pieces to be welded with plaster or low concentration covering. As there is not heating it is not necessary to use high temperature covering.
- Also use, for temporary fixing, pieces or bars of the base metal which will be removed after the welding.
- It is usually necessary to add some metal for bonding bridges after cutting. Use carbon-free wire from the same base metal with a diameter of 0.5mm to 0.7mm. Ask your metal supplier for it.
- Chamfer the edges and put the adding metal between the piece and the electrode. Repeat the operation until the complete bonding of the parts. To finish, double the pulse and depth adjustment.
- Repairing holes and porosity. Weld it with regulation point at 3 on the pulse and also 3 on the depth over the damaged area to open it. Next, repair the defect with adding metal.
- When welding it is possible that there is some distortion on the welding area. For repairing, weld also on the opposite direction.
- With dental braces it is not usually necessary to add metal.

### 4.8 Tungsten electrode

For a good weld result, it is necessary that the electrode is always sharpened with 25' angle. It is common that its tip gets round on the first welds due to inexperience. Use fine-grain carborundum or the diamond disc to sharpen the electrode. We usually sharpen both tips.


**Warning:** Use only the electrodes supplied by EDG. Electrodes made with different material or of bad quality will severely damage the equipment.

### 4.9 Grounding wires

To form the plasma, it is necessary that the pieces to be welded are grounded, to do so use the grounding wires with clips (2) Pic. 6.

- Connect the wires to the grounding terminals (1) Pic. 6.
- Connect the wire clip to the piece. One clip is usually enough to ensure good contact. When the pieces which will be welded do not have contact it is necessary to use one clip for each section.

## 4.10 Cooling down

When welding large structures or the ones with resins or ceramic, to keep the temperature under control we can cool them down by placing the piece in front of an air nozzle (3) Pic. 6 and activate the switch (4) Pic.6. to do so it is necessary that the compressed air is connected (item 4.4).

## 4.11 Hand rests


The hand rests are very useful to keep the hands stable (5) Pic. 6.

## 4.12 Control panel

- 1- Main switch. Turns the equipment on and off.
- 2- Contact buzzer key. Exclusive system which sounds an alarm when the piece to be welded touches the electrode making the operator's job easier.
- 3- Pulse time adjustment. It adjusts the time in which the plasma is produced.
- 4- Pulse time digital display (milliseconds)
- 5- Penetration depth adjustment. It adjusts the alternate current intensity that will form the plasma and consequently the welding depth.
- 6- Penetration depth digital display (related to amperage).

## 4.13 Magnifying display window, photo sensitive crystal

Micromelt is equipped with tilting lens with 3x zoom and photo sensitive crystal that activates automatically a dark filter when the plasma is formed protecting the operator's eyes.

**Warning:** The plasma generates ultraviolet light which is harmful for the eyes, so never use the equipment without the photo sensitive filter or with the flap door open.

## 5 Operation

### 5.1 Argon outflow adjustment

- The best argon outflow ranges from 3 to 5 liters/min (0.1 to 0.2ft<sup>3</sup>/min), to adjust it:
- Turn the main switch on.
  - Turn the argon pressure regulator knob all the way to the left until there is no mere resistance.
  - Open the argon cylinder valve.
  - Keep the pedal activated and turn the pressure regulator knob until the outflow indicator gauge shows 5 liters/min (0,2 ft<sup>3</sup>/min).
  - Release the pedal, the argon will keep flowing for a few seconds.
  - Check if there are argon leakages on the hose connections using soap foam.

### 5.2 Welding

- Place the pedal in a comfortable position.
- Open the argon cylinder valve.
- Turn the main switch on.
- Turn the buzzer key on.
- Turn the pulse button until the display show "05".
- Turn the depth button until the display shows "03".
- Fix to one of the clips a chrome nickel or chrome cobalt structure which has not been used.
- Use the hand rests for more stability.
- Push the pedal. You will hear the argon flowing.

- Put the structure against the electrode without pressing keeping the hands firm, the buzzer will indicate the contact. Then release the pedal, the weld will be done