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
--

1 Receiving

The equipment is packed separately from the vacuum pump, check if you have two boxes if you have acquired the furnace and the pump. Check the general conditions of the package, if any damage is noticed contact the shipping company. We remind you that the equipment is shipped at the purchaser's account and risk and it is insured by the shipping company.

notations

1.1 The Sinter Press package should include

- A) A furnace settled in injected foam
- B) A spool for conventional ceramic Pic4-Item-15;
- C) A spool for pressed ceramic. Pic-5-item-19.
- D) A rigid mat -I for ceramic. Pic-4-Item-17.
- E) A rigid mat -II for sintering Pic-4-Item-16.
- F) A soft mat.
- G) A refractory platform for pressed ceramic Pic-5-Item-20.
- H) A see through refractory platform for ceramic.
- I) A set of isothermal pins. Pic-4-Item-18.
- J) Three 1.5m hoses N. A. T.
- K) Metallic clamps. (6 pieces).
- L) Instruction manual.
- M) Warranty document.

1.2 The vacuum pump package should include

- A) The vacuum pump
- B) Rubber, floor insulating caps for fixation (4 pieces).
- C) Connectors identification labels

2 S.A.L.V. System

Your Sinter Press has highly technological devices that allow it to perform 4 functions perfectly: traditional ceramic burn, pressed ceramic, sintering and infiltration of alumina composites. When sintering alumina, the lining used heated over 500° C releases sub-products which are absorbed by the muffle under vacuum and can contaminate the ceramic, staining and fissuring it during the sintering process.

Long periods under high temperatures (1,200°C) in pressurized chambers (vacuum muffles) can shorten the heating resistance useful life. The S.A.L.V. system was specially developed to overcome such problems

Lift					Approximate Time
S	M	.	D	M	Sai do fim de curso
S	D	.	D	D	Percurso contínuo
S	1	.	D	1	1 Minute
S	2	.	D	2	2 Minutes
S	3	.	D	3	3 Minutes
S	4	.	D	4	4 Minutes
S	5	.	D	5	5 Minutes
S	6	.	D	6	6 Minutes
S	7	.	D	7	7 Minutes
S	8	.	D	8	8 Minutes
S	9	.	D	9	9 Minutes

17-18) KEYS (.) PONTO, 0,1,2,3,4,5,6,7,8,9 :- They allow the insertion of new data for the recipes variables. Activate the desired key and insert the new value.

17-19) KEYS WITH SYMBOLS:- They allow the insertion of customized identification for the recipe names during the editing and data saving stages.

Specification:

- State-of-the-art micro processed control MOTOROLA.
- SMD electronic circuit cards
- Liquid crystal display with BACKLIGHT. Blue background
- Power = 1800 Watts.
- Power supply: 110 or 220 volts. 50/60 Hz. Under request
- Fuses: 20 amperes to 110 volts, 15 amperes to 220 volts.
- RS-232 communication – for technical assistance.
- Alphanumeric keyboard with 39 programmes and monitoring keys.
- Maximum operation temperature: 1200°C
- Nominal working pressure 60 Lbf / in² (Piston)
- Nominal vacuum work -27 in Hg. / mm Hg.

Furnace dimensions:

- Width.....330mm
- Height.....800mm
- Depth.....420mm
- Net weight..... 22.5Kg
- Gross weight.....28.5Kg

Package dimensions:

- 440mm x 840mm x 490mm

SALV®

Sistema de aquecimento de longa vida.
(Long Life Heating System)

Two years guarantee at high temperatures.

The metallic alloys used as heating resistance have chemical elements in its composition that, when heated and in contact with air, create an oxide adherent layer which protects them from outwear.

In ceramic furnace muffles, the process for the oxide layer formation can be impaired due to the resistance imprisonment and to the vacuum operation.

The exclusive S.A.L.V. system, developed and patented by EDG, promotes airflow inside the quartz pipe that supports the muffle resistance. This flow renews and maintains the protective oxide layer, increasing the useful life of the resistance at high temperatures.

This exclusive technological advance allowed EDG to extend the guarantee to 2 years for the furnaces equipped with this system.

To endure the high temperatures involved in the processes, the sensor (termopar) should be made of platinum which, on the other hand, is not adequate for low temperatures and for the necessary speed for eliminating water during the sintering process. High level software was developed to compensate the platinum limitations in low temperatures. Due to the characteristics of the sintering process, new thermal insulating materials were used and an efficient ventilation system keeps the components temperature under adequate values.

3 Main Characteristics

- Quartz muffle with thermal insulation ceramic mat molded through low thermal mass vacuum, which insures perfect temperature homogeneity.
- The exclusive system S.A.L.V® Sistema de aquecimento de longa vida (long life heating system)
- Muffle forced ventilation for the sintering/infiltration process which eliminates stains and minimizes fissures.
- 99 programs: 79 for ceramic and 20 for sintering/infiltration.
- Management by microprocessing.

- Easy operation and visualization of the burning parameters.
- Ambient working temperature at 1,100°C for ceramic and at 1,200°C for sintering/infiltration.
- Maximum temperature limiter at 1,200°C
- Linear decimal heating speed from 0.1 to 70.0°C/min.
- Thyristored power control
- Ceramic burn time from 0 to 30 minutes.
- Sintering/infiltration time from 0 to 9 hours 59 minutes.
- Decreasing time visualization.
- Vacuum time from 0 to 30 minutes with decreasing visualization.
- Automatic or manual (which can be activated ant any stage of the process) vacuum operation
- Programmable and independent ascending and descending elevator times.
- Safety systems protect the equipment from operation mistakes or defects.
- Flashing and beeping signs of all process steps.
- To avoid ambience heating and save energy, an alarm will sound after 5 minutes if the equipment is not being operated the muffle will be partially closed automatically.
- Night Mode: after 2 hours without any operation, the muffle is closed and the temperature is kept at 100°C.
- Programmable quick cooling.
- Operation time totalizer.
- Low cost muffle replacing system.
- Elevator with bottom entrance perfectly balanced
- Electronic system forced ventilation keeps the components safe.
- ARGON SYSTEM, allows the cleaning and helps with ceramic burn emphasizing its colors, enhancing the translucence and opalescence.
- Serial communication for maintenance and system checking.
- Temperature safety devices, temopar, vacuum, elevator, power, muffle with S.A.L.V system, piston and compressed air system.
- Digital panel in liquid crystal with 4 lines by 20 columns on a blue background.
- Sensitive membrane with sound recognition of the operator's touch.

4 Installing

Note: The numbers in parentheses refer to the picture that follows.

17-11) KEY - START STOP:- This key has two functions. When it is activated with the muffle open it will perform the burn START; when it is activated with the muffle closed it will perform the STOP, interrupting the process.

17-12) KEY - ARGON:- This function allows the automatic argon washing at the beginning of the burn process. Activate the PROGRAMA key until you reach (tela 4) and confirm the function ARGON with the SETA ESQUERDA (left arrow key), filling up the respective field.

17-13) KEY - RESFR:- This function helps forced cooling to occur between burn processes using the vacuum pump to make cold air circulate in the muffle helping the heat exchange. Activate the PROGRAMA key until you reach (tela 4) and confirm the function RESFR with the SETA ESQUERDA (left arrow key), filling up the respective field.

17-14) KEY - NIGHT:- This function allows the furnace to be resting in economic mode, while it is not being used. Ex. during the night it keeps the muffle at 100°C.

17-15) KEY - AUTO START:- This function allows performing a new burn process without the need to reach T1. Ex. After turning the furnace on, program it and place the piece to burn and activate AUTO START. The process will be automatically performed.

17-16) KEY - SOBE:- This function allows the elevator electric activation using the SOBE (ascent) key at the operator's touch. When it is activated it will be on until the next touch to interrupt the operation. It is turned off automatically when it reaches the end of high course. To program drying values use the PROGRAMA (tela-5) key and type the desired value. Check the values table below.

17-17) KEY - DESCE:- This function allows the elevator electric activation using the DESCE (descent) key at the operator's touch. When it is activated it will be on until the next touch to interrupt the operation. It is turned off automatically when it reaches the end of low course. To program piece slow exit use the PROGRAMA (tela-5) key and type the desired value. Check the times table below.

17-2) KEY - PROGRAMA:- Allows you to program and navigate through the necessary screens for the burn. Each touch on this key will be visualized on the display, type all desired information and activate the PROGRAMA key again to go to the next programming stage, 5 sequential screens should be filled in to complete the procedure.

17-3) KEY –MONITORAÇÃO:- This key allows you to visualize the burn variables during the process on a single screen in which the data will be updated while the stages are performed.

17-4) KEY - A/N:- This key allows you to alter the keyboard writing mode between numeric and alphabetic according to the programming needs when editing the recipes names.

17-5) KEY – LIMPA REC :- This key can be used when editing and programming recipes to erase all content on a screen. EX- T2, PAT2, VAC2, VEL2. or when we want to exclude a recipe from the memory. This function is only possible when in GRAVAÇÃO (saving) mode.

17-6) KEY – LER:- Activate the LER (read) key to load a recipe from the control memory, to do so type the desired number and activate the LER (read) key again. The operation will be then concluded.

17-7) KEY - GRAVA:- Activate the GRAVA (save) key to save a new recipe in the control memory. Choose the number, the name and activate the GRAVA (save) to conclude the operation.

17-8) SETAS (arrow keys):- Use the direita (right) and esquerda (left) arrow keys to navigate through the selected digits during programming stage, and also to activate the argon, cooling and auto-start fields filling in.

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
17-10) KEY - ENVIA : This key is used when we want to send a new recipe to the furnace manually or when only one datum is altered. It can be sent to be performed by the furnace. We remind you that the ENVIA (send) function is automatically performed every time we go to the monitoring screen.

Your Sinter Press should be far from curtains and inflammable material. A furnace is a heat generator which needs to be dissipated; otherwise its components will be overheated. Therefore, place the equipment at a ventilated area which allows free air circulation.

It is advisable to leave a 15cm clearance between the furnace and anything that can impair the ventilation.

Place the furnace far from taps or basins which can sprinkle water over the equipment.

Pic. 1

4-1) Check if your network voltage is the same as the one shown on the furnace Pic (1).- item (1).

4-2) Install the furnace into an exclusive electrical network using 6mm wiring to a 110 volts network, or 4mm wiring to a 220 volts network.

4-3) Never connect the furnace to the same network in which other furnaces, compressors, electrical taps or any other high electricity consumption devices are already connected.

4-4) Check if the outlet in which the furnace will be plugged in is in perfect condition and of good quality to bear at least 20 amperes.

4-5) Connect the grounding end (round pin) to a grounding rod and never to the electrical network neutral.

Note: Not checking the procedures above will interfere on the equipment correct working, on its guarantee and also on the operator safety.

4-6) If your electrical network is very instable (the lights flash a lot), it is advisable to use a voltage regulator such as a saturated autotransformer with at least 2kw capacity.

Warning: Never use, by any means, regulators used in computers.

4-7) Connect the vacuum hose to the furnace Pic-1 item (5). Plug the vacuum pump power cable to the outlet which is at the backside of the furnace Pic-1 item (7).

4-8) connect the compressed air hose to the furnace Pic-1 item (11). The air should be supplied dry by the compressor at a minimum pressure of 2 kg/cm² (30lbf/pol²). This pressure should be kept while the furnace is working with pressed ceramic.

Warning: The lack of pressure during this process will abort the burn and a message will be shown on the control display.

4-9) Connect the argon hose to the furnace and to the cylinder outflow regulator Pic-1 item (6). The argon used is of commercial type. The cylinder and the regulator are not supplied with the equipment. We recommend to contact a gas supplier in your area. The maximum argon outflow should be adjusted in 10 L/m (liters per minute). Note: We remind you that the argon provides better colors, opalescence and translucence. However, its use can be discharged. For that do not connect the argon hose and do not program the function for argon. Check the hoses connections if they are tight to avoid unpleasant leakages.

Message:- Temperature over the allowed limit – check – (OK)

Error - 11 - Piston flaw ** COMPRESSOR **

Message:- Low compressed air pressure – check – (OK)

Error- 12 - Programming flaw ** PRENSAGEM **

Message:- Pressing time should be lower than the burn time – check – (OK)

Error - 13 - Programming flaw ** PRENSAGEM **

Message:- Vacuum time and burn time should be the same – check – (OK)

Error - 14 - Programming flaw ** ARGÔNIO **

Message:- It is not possible to program the argon washing without vacuum – check – (OK)

Error - 15 - Programming flaw ** NIGHT **

Message:- This function can only occur in pat. T1 with the elevator on the platform – check – (OK)

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Message:- It was not possible to find the furnace platform. – check – (OK)

Error - 17 - Movable tray device flaw ** START **

Message:- Elevator tray out of the correct platform – check – (OK)


Error - 18 - Technical access flaw ** ROT-TÉCNICA **

Message:- Access granted only for authorized personnel – check – (OK)

Error - 19 - Operator access flaw ** SENHA **

Message:- password does not match, try again – check – (OK)

17 Keyboard Functions

17-1) The alphanumeric keyboard, has 39 keys for functions programming and burn system monitoring which are divided in blocks, numeric, alphabetic, auxiliary functions and control.

being saved or sent to the furnace. When this happens, wrong information is not carried on since they could harm the burn process.

15-4) System flaw messages: electric problems may occur and interfere with the equipment best performance. To avoid flaws some vital functions are continually monitored and if there is any error, it will be signalized.

16 Messages definitions

16-3) ERRORS MODE: - They are messages which signal general problems; remain on the screen until the operator fixes the problem, releasing the operation then with the OK key.

Error 01- thermocouple burst-out flaw ** TERMOPAR **

Message: Thermocouple was damaged, contact technical assistance – check – (OK) The process being performed is aborted, an alarm will sound warning about the flaw. Contact the technical assistance.

Error 02 – Elevator flaw. ** ELEV-SOBE **

Message: electromechanical problems – check – (OK)

Error 03 - Elevator flaw ** ELEV-DESCE **

Message: electromechanical problems – check – (OK)

Error 04- Velocity flaw ** VELOCIDADE **

Message:- It is not possible to save a ZERO value for T2 velocity – check – (OK)

Error 05- Velocity flaw ** VELOCIDADE **

Message:- It is not possible to save a ZERO value for T1 velocity – check – (OK)

Error - 06 – Saving flaw - ** REC-PADRÃO **

Standard recipes cannot be overwritten – check – (OK)

Error - 07 – Vacuum flaw – ** FALHA-VÁCUO **

Mensage:- It was not possible to produce the vacuum – check – (OK)

Error - 08 – Vacuum flaw – ** FALHA-VÁCUO **

Mensage:- It was not possible to turn the vacuum off – check – (OK)

Error - 09 - Version updating flaw ** PROGRAMA **

Mensage:- Version incompatible with the current furnace platform – check – (OK)


Error - 10 - A/D conversion flaw ** TEMPERATURA **

5 Control panel

5-1) The panel has liquid crystal display with 4 lines by 20 columns on a blue background and yellow printing with BACKLIGHT inside illumination that allows easy identification of the programming functions and log monitoring.

5-2) It has a membrane sensitive to the operator's touch, which makes the functions activation easier with beeping sound recognition.

5-3) Panel identification and its main functions. PIC-2

5-4) VISOR (display) – Liquid crystal display 4 / 20 blue. Item-21.

5-5) PROGRAMA (program) – Sequential program key. Item-22.

5-6) MONITORAÇÃO (monitoring)– control system monitoring key. Item-23.

5-7) FUNÇÕES RÁPIDAS (quick functions) – Use the PROGRAMA key, the sequential screens for programming and monitoring - Item-24.

5-8) ELEVADOR SOBE E DESCE (elevator up and down)– touch sensitive digital command which turns the elevator motor on and off. Item-25.

5-9) FUNÇÕES AUXILIARES (auxiliary functions) – Argon washing, quick cooling, pre burn and pressing. Item-26.

5-10) NIGHT – Automatic or manual resting system, it only operates with ceramic. Item-27.

5-11) AUTO-START – This operation allows the automatic burn start from the ambient temperature. Item-28.

5-12) START/STOP - It allows to start or interrupt the burn through the operator's

command. Item-29.

5-13) CONJUNTO DE TECLAS COM FUNDO CINZA (Keys set with gray background)– Programming and control system exclusive operation. Item-30.

5-14) TECLADO NUMÉRICO (numeric keyboard) – Keyboard for inserting the numeric content for the recipes variables. Item-31.

5-15) MENSAGENS (messages) – Field reserved for visualizing the process stages. Item-32.

6 System start-up – step-by-step

6-1) Plug the furnace to a good quality wall socket properly grounded to a copper rod and never to the network neutral, check if your network voltage is the same as the one shown on the furnace label. Pic-1 item (1).

6-2) Turn the main switch on Pic-1 item (3) and wait for the equipment start-up. The first screen on contains the program versions data. Ex-1.

STATUS			REC°N		TEMP °C			VÁCUO		
E	D	G	E	Q	I	P	A	M	E	N
E	D	G	H	A	N	D	V	2	.	2
E	D	G	B	A	S	E	V	1	.	3
W	W	W	.	E	D	G	.	C	O	M
ELEVADOR			MENSAGENS							

6-3) The next screen shows the starting-up and the whole system checking-up, a scroll bar keeps running through the screen while it performs the initial configurations. Ex-2.

STATUS			REC°N		TEMP °C			VÁCUO		
I	N	I	C	I	A	L	I	Z	A	N
I	N	I	C	I	A	L	I	Z	A	N
I	N	I	C	I	A	L	I	Z	A	N
I	N	I	C	I	A	L	I	Z	A	N
ELEVADOR			MENSAGENS							

6-4) When the checking-up is over, the system loads the sintering standard recipe nº 80 and waits for the operator's command to perform a new function. We remind you that when sintering, the furnace will just heat up after the start command is activated; therefore the temperature on screen is close to the ambient temperature. Ex-3.

STATUS			REC°N		TEMP °C			VÁCUO		
M	O	N	I	T	0	8	0	0	0	2
T	1	=	0	2	0	0	P	A	T	1
V	E	L	1	=	0	5	.	H	0	:
S	D	D	M	.	P	A	D	R	Ã	O
ELEVADOR			MENSAGENS							

to protect the equipment, avoiding energy waste and keeping the muffle heated to a minimum temperature of 100°C, avoiding humidity. Note: To deactivate this function activate the NIGHT key again.

Note: To deactivate this function activate the NIGHT key again.

14-5) AUTO START:- This option allows the automatic burn start from the moment the temperature reaches the T1 baseline with the elevator on the platform. To do so it is necessary that the piece is already placed on the burn tray. A beeping and visual alarm will signal the process start.

Note: The burn will be performed once.

This option can be activated either on the T1 ascent rampa from 25° C to 600 ° C or on the descent from 960° C to 600° C when cooling, activating the AUTO-START key, or it can also be programmed through the PROGRAMA key - (tela 4), in which it will make part of the recipe. EX-55.

Programming auxiliry functions.

STATUS			REC°N		TEMP °C			VÁCUO		
T	E	L	A	4	0	6	5	A	R	G
P	R	É	Q	=	0	5	:	R	E	S
P	R	E	N	=	0	0	:	S	F	R
N	O	M	E	:	.	.	.	A	U	T
ELEVADOR			MENSAGENS							

PROGRAMANDO
PRÉ-QUEIMA
ARGON
RESFRIAMENTO
E AUTO-START
EX - 55

15 Alarms and messages

15-1) Most messages shown on the display aim to indicate some stages during the process operation, however programming errors indication are also very common. There are also protection messages for the control system which will accuse electronic circuits operation flaws.

15-2) Operation messages: - All functions which are being performed are shown on the display on a specific field. Each one of them is shown according to the burn current stage, for example: - RAMPA (ramp), PATAMAR (baseline), PRÉ QUEIMA (pre burn), RESFRIAMENTO (cooling), AUTO START, among others, these messages aim to illustrate the burn process during it is being performed.

15-3) Programming messages:- They accuse flaws on the programming which is

STATUS		REC°N		TEMP °C		VÁCUO											
E	R	R	O	=	0	7	*	*	*	*	V	Á	C	U	O	*	*
N	Ã	O	F	O	I	P	O	S	S	I	V	E	L				
F	O	R	M	A	R	V	Á	C	U	O	.						
						.V	E	R	I	F	I	Q	U	E	(O	
ELEVADOR		MENSAGENS															

- Start flaw. EX-93 – ERRO=17

STATUS		REC°N		TEMP °C		VÁCUO											
E	R	R	O	=	1	7	*	*	*	*	S	T	A	R	T	*	*
B	A	N	D	E	J	A	D	O	E	L	E	V	A	D	O		
F	O	R	A	D	A	B	A	S	E	C	O	R	R	E	T	A	
						.V	E	R	I	F	I	Q	U	E	(O	
ELEVADOR		MENSAGENS															

13 Quick functions

These functions allow the programming or monitoring of the desired variable through the PROGRAMA key shortcut.

14 Auxiliary functions


14-1) AUXILIARY FUNCTIONS:- They are complementary resources that help the burn process, improving the final work quality. They make the control system operation easier, protecting the equipment and controlling energy waste, optimizing total process time. The auxiliary functions are: ARGON SYSTEM, QUICK COOLING, AUTO-START, PRE BURN, NIGHT, EX-55.

14-2) ARGON SYSTEM:- This function performs the automatic piece washing at the beginning of the burn process. Use only inert and non oxidant argon gas, it helps emphasizing the colors, enhancing translucence and opalescence. EX-55.


14-3) QUICK COOLING:- This function is used for optimizing the time spent between burn processes, forcing the heat exchange to T1 allowing a new burn process to start. EX-55.

14-4) NIGHT:- This function is either automatically activated when the equipment is not being used for a two-hour period or manually by using the NIGHT key, which puts the furnace under complete resting, retreats the elevator and lowers the T1 temperature to 100 ° C. This option can only be activated with T1 baseline for ceramic burn with the elevator on the platform. The main objective of this function is


6-5) The furnace should not start the burn process without the spool on the tray; otherwise the sealing ring and the tray will be seriously damaged and will need to be replaced. Choose the proper accessory for the type of work to be burn. PIC-4 e PIC-5.

Pic.4

Pic.5

Pic.6

6-6) Accessories for the burning - SINTERING.

- A) Use spool - item (15).
- B) Use the rigid mat item (16) on the spool.

6-7) Accessories for the burning - .CERAMIC – Conventional.

- C) Use spool – item (15)
- D) Use the rigid mat item (17) on the spool.
- E) Use isothermal pins for supporting the piece. Item (18).

6-8) Accessories for the burning -.CERAMIC – PRESSED.

- F) Use spool – item (19).
- G) Use refractory platform item (20) inserted in the spool to place the ring.

Pic.3

Warning: The set of RINGS, PADS and PLUNGER should not be higher than the maximum height of 92 millimeters, as shown on PIC-(6).

6-9) Move the elevator tray to the front item-(14) holding it by its knob and place the spool item-(15) on its fitting, next bring the tray back to its correct position on the platform to start the equipment operation.


7 Programming ceramic recipes

7-1) For better illustrating, we will follow some practical recipes examples for conventional ceramic, pressed ceramic, sintering and infiltration.

Ex-1 recipe model- conventional ceramic

EXEMPLO DE RECEITA PARA CERÂMICA CONVENCIONAL										
TEMP-1 °C	VEL-1 Fixo	PAT-1 Elevador	TEMP-2 °C	VEL-2 °C/min	PAT-2 min	VAC-2 min	TEMP-3 °C	VEL-3 °C/min	PAT-3 min	VAC-3 min
600	60	BASE	960	55	01:00	01:00	980	20	01:00	0

PRÉ. QUEIMA min	PRENSAGEM min	ARGON		RESFRIAMENTO		AUTO-START Campo preenchido	ELEVATOR	
		Campo preenchido	Campo preenchido	Campo preenchido	Campo preenchido		UP	DOWN
05:00	00:00	X	YES		NO		5	3

7-2) Start the programming by reading the standard ceramic recipe (01); it will serve as reference for programming the other recipes.

Process messages

MENSAGENS									
.	L	A	V	A	G	E	M	A	R
.	.	R	E	S	F	R	I	A	G
.	N	I	G	T	O
I	N	I	C	I	A	R	Q	U	E
.	A	M	P	A	M
.	P	A	T	A	A
.	R	E	N	S	R
.	P	R	E	N	S
.	R	E	N	S	R
.	A	U	T	O	S
R	E	S	F	.	A	U	T	O	S
.	S	T	O	P	.
.	S	O	B	E	.
.	D	E	S	C	E
.	P	R	O	T	E
.	A	D	R	Ã	O
S	U	B	S	T	P	A	R	Ã	O
N	O	M	E	S	R	E	C	E	I
.	R	E	C	E	I
N	O	M	E	S	E	C	E	I	T
.	T	A	(S/N)	.	A
N	O	M	E	S	R	E	C	E	I
.	E	I	T	A	S

12-2) Variables such as T1, T2, T3, Vel1, Vel2, Vel3, Pat1, Pat2, Pat3, Vac2, Vac3, Pre burn, Pressing, recipe number, elevator ascent(D) e descent(D), are continually updated during the burn process on the monitoring screen. EX-38.

Note: The saved recipes names will be shown on the first stage of the burn process, ambient temperature ramp until it reaches the T1 baseline.

Emphasizing items are on this screen EX-38

STATUS	REC°N	TEMP °C	VÁCUO	MONITORAÇÃO CERÂMICA EX - 38													
				M	O	N	I	T									
M	0	0	2	0	6	0	0	-	2	7							
T	2	=	0	9	6	0	P	A	T	2	=	0	3	:	0	0	
V	E	L	2	=	0	4	0	V	A	C	2	=	0	2	:	0	0
S	D	D	D	.	.	.	R	A	M	P	A
ELEVADOR		MENSAGENS				MENSAGENS				MENSAGENS				MENSAGENS			

12-3) System flaws monitoring. This piece of equipment has several safety devices which permit optimal working conditions, if there is any kind of flaw the burn is aborted and an error message will remain on the screen until the operator fixes the problem, releasing the operation then with the OK key. Some examples for the most common errors follow:

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11-6) On the next screen choose recipe 81 typing number 81 on the REC field or navigate through the options list using the left or right arrows, the recipe name will be on the message field [. . SINTERIZAÇÃO . .] if you want to overwrite a recipe, confirm the operation activating the GRAVA (save), in this case a message will be shown [SUBST. RECEITA? (S/N)] if you want to overwrite, choose SIM (yes), type the letter [S] and the operation will be concluded and the recipe will be sent to the platform. If you choose not to overwrite, choose NÃO (no) and type the letter [N]. The cursor will return to the REC=002 field for changing the recipe number. To change the recipe name, activate the OK key and the cursor will go to the field reserved for the recipe name, type the new name, EX. (. . INFILTRAÇÃO . .) confirm the operation again using the GRAVA (save) key. The recipe will be automatically sent to the furnace platform. EX-37.

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11-7) After reading the recipe the system will return to the monitoring screen. EX-33.

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12 System monitoring

12-1) Burn process monitoring. This screen allows visualizing all burn variables during its process. Each stage of the burn process is shown on the display message field, making the operation monitoring easier.

7-3) Activate the LER (read) key and choose the type of work to be performed, ceramic or sintering. To alternate between the two options use the OK key and confirm the operation again using the LER (read) key. Ex4.

STATUS		REC°N	TEMP °C	VÁCUO	
L	E	R			
0	1	à	7	9	
8	0	à	9	9	
		C E R	Â M I C A	S I N T E R I Z A Ç Ã O	
					LEITURA DAS RECEITAS DE CERÂMICA
					EX - 4
		ELEVADOR	MENSAGENS		

7-4) On the next screen choose receita 01 typing número 1 on the REC field, or navigate through the options list using the esquerda (left) or direita (right) arrows, the name of the recipe will be shown on the messages field (PADRÃO CERAM) EX-5, at this time you should confirm the option by activating the LER (read) key again to complete with the operation. The recipe will automatically be sent to the platform and when there is no recipe recorded, the message (RECEITA VAZIA) will be on the display.

STATUS		REC°N	TEMP °C	VÁCUO	
L	E	R			
R	E	C :	0	0	1 . . C E R Â M I C A . .
A	L	F	A		
N	O	M	E :	.	P A D R Ã O . C E R Â M .
ELEVADOR		MENSAGENS			LEITURA DAS RECEITAS DE CERÂMICA
					EX - 5
		ELEVADOR	MENSAGENS		

7-5) After the standard recipe is read the system will return to its monitoring screen. EX- 11

STATUS		REC°N	TEMP °C	VÁCUO	
M	O	N	I	T	
M	O	N	I	T	0 0 1 0 6 0 0 0 0 0
T	E	L	I	1	= 0 6 0 0
S	D	D	M	1	= 0 6 0 P A D R Ã O . C E R Â M .
ELEVADOR		MENSAGENS			MONITORAÇÃO
					EX - 11
		ELEVADOR	MENSAGENS		

7-6) To start programming a new recipe you should activate the PROGRAMA key, and the display will show the message TELA1, the cursor will be flashing over the field related to temperature T1=600. At this time we have two options, we can simply type new values on the already filled fields or we erase all values using the LIMPA REC key, then we can type on the empty fields. EX-12

Note: For CONVENTIONAL CERAMIC EX-1, we will use the empty fields activating the function LIMPA REC.

STATUS		REC°N	TEMP °C	VÁCUO	
T	E	L	A	1	0 0 1
T	1	=	0	0	0 0 0
N	O	M	E	:	.
ELEVADOR					MENSAGENS

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EX - 12

7-7) On the screen TELA-1, now with the empty fields, you will see the cursor flashing on the T1=0000 field. Type the new value T1=0600 following the steps presented for standard conventional ceramic recipe EX-6.

STATUS		REC°N	TEMP °C	VÁCUO	
T	E	L	A	1	0 0 1
T	1	=	0	6	0 0
N	O	M	E	:	P A D R Ã O . C E R Ã M
ELEVADOR					MENSAGENS

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EX - 6

7-8) Activate the PROGRAMA key again and a new screen TELA-2 can be visualized, in which the cursor will be flashing on the T2=000 field. Type the new value of T2=960 and confirm using the OK key, then the cursor will be flashing on the velocity field VEL2=000, type the new value of VEL2=55 and confirm it using the OK key. Now type the new value for the baseline time PAT2=01:00 and vacuum time VAC2=01:00, these times are programmed using minutes. Confirm the operation by using the OK key. EX-7.

STATUS		REC°N	TEMP °C	VÁCUO	
T	E	L	A	2	0 0 1
T	2	=	0	9	6 0
V	E	L	2	=	0 5 5
N	O	M	E	:	P A D R Ã O . C E R Ã M
ELEVADOR					MENSAGENS

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EX - 7

7-9) Activate the PROGRAMA key again and a new screen TELA-3 can be visualized, in which the cursor will be flashing on the T3=000 field. Type the new value of T3=980 and confirm using the OK key, then the cursor will be flashing on the velocity field VEL3=000, type the new value of VEL3=020 and confirm it using the OK key. Now type the new value for the baseline time PAT3=01:00 and vacuum time VAC3=00:00, these times are programmed using minutes. Confirm the operation by using the OK key. EX-8

STATUS		REC°N	TEMP °C	VÁCUO	
T	E	L	A	3	0 0 1
T	3	=	0	9	8 0
V	E	L	3	=	0 2 0
N	O	M	E	:	P A D R Ã O . C E R Ã M
ELEVADOR					MENSAGENS

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EX - 8

11-2) EX-35 – Opaque burn. Activate the GRAVA (save) key and choose the type of work to be performed, ceramic (cerâmica) or sintering (sinterização). To alternate between the two options use the OK key and confirm the operation again using the GRAVA (save) key. Ex.34

STATUS		REC°N	TEMP °C	VÁCUO	
G	R	A	V	A	C E R Â M I C A
0	1	=	7	9	S I N T E R I Z A Ç Ã O
8	0	=	9	9	
N	O	M	E	:	ELEVADOR
					MENSAGENS

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EX - 34

11-3) On the next screen choose recipe 02 typing number 02 on the REC field or navigate through the options list using the left or right arrows, the recipe name will be on the message field [. . . OPACO . . .] if you want to overwrite a recipe, confirm the operation activating the GRAVA (save), in this case a message will be shown [SUBST. RECEITA? (S/N)] if you want to overwrite, choose SIM (yes), type the letter [S] and the operation will be concluded and the recipe will be sent to the platform. If you choose not to overwrite, choose NÃO (no) and type the letter [N]. The cursor will return to the REC=002 field for changing the recipe number. To change the recipe name, activate the OK key and the cursor will go to the field reserved for the recipe name, type the new name, EX. (. . . OPACO . . .) confirm the operation again using the GRAVA (save) key. The recipe will be automatically sent to the furnace platform. EX-35.

Saving opaque burn. REC-02.

STATUS		REC°N	TEMP °C	VÁCUO	
G	R	A	V	A	C E R Â M I C A
R	E	C	: 0	0	2
A	L	F	A		
N	O	M	E	:	O P A C O
ELEVADOR					MENSAGENS

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EX - 35

11-4) After reading the recipe the system will return to the monitoring screen. EX-30.

STATUS		REC°N	TEMP °C	VÁCUO	
M	O	N	I	T	0 0 2
T	1	=	0	6	0 0
V	E	L	1	=	0 6 0
S	D	D	D	=	O P A C O
ELEVADOR					MENSAGENS

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EX - 30

11-5) EX-37 – Sintering burn. Activate the GRAVA (save) key and choose the type of work to be performed, ceramic (cerâmica) or sintering (sinterização). To alternate between the two options use the OK key and confirm the operation again using the GRAVA (save) key. Ex.36

STATUS		REC°N	TEMP °C	VÁCUO	
M	O	N	I	T	0 0 3
T	1	=	0	8	0 0
V	E	L	1	=	0 6 0
S	D	D	D	.	P R E N S A D A . . .
ELEVADOR		MENSAGENS			

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10-8) EX-32 – Sintering burn. Activate the LER (read) key and choose the type of work to be performed, ceramic (cerâmica) or sintering (sinterização). To alternate between the two options use the OK key and confirm the operation again using the LER (read) key. Ex.19

STATUS		REC°N	TEMP °C	VÁCUO	
L	E	R			
0	1	à	7	9	
8	0	à	9	9	
			C E R	Â M I C A	
			S I N T E R	R I Z A Ç Ã O	
ELEVADOR		MENSAGENS			

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10-9) On the next screen choose recipe 81 typing number 81 on the REC field or navigate through the options list using the left or right arrows, the recipe name will be on the message (. . . INFILTRAÇÃO . . .) At this stage, confirm the operation activating the LER key again to complete it. The recipe will be automatically sent to the furnace platform. EX-32.

Sintering and infiltration burn. REC-81

STATUS		REC°N	TEMP °C	VÁCUO	
L	E	R			
R	E	C :	0 8 1	S I N T E R I Z A Ç Ã O	
N	O	M	E :	. . . I N F I L T R A Ç Ã O . . .	
ELEVADOR		MENSAGENS			

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10-10) After reading the recipe the system will return to the monitoring screen. EX-33. Sintering monitoring REC-81

STATUS		REC°N	TEMP °C	VÁCUO	
M	O	N	I	T	0 8 1
T	1	=	0	2	0 0
V	E	L	1	=	0 5 . 0
S	D	D	M	.	I N F I L T R A Ç Ã O . . .
ELEVADOR		MENSAGENS			

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11 Saving recipes

11-1) Saving recipes in the database.

Ex-35 - Conventional ceramic opaque burn. REC-02

Ex-37 - Sintering and infiltration burn. REC-81

7-10) Activate the PROGRAMA key again and a new screen TELA-4 will be visualized, in which the cursor will be flashing on the pre burn field PRÉQ=00:00, type the new value of PRÉQ=05:00 and confirm it using the OK key, then the cursor will be flashing on the press field PREN=00:00, keep the value in ZERO and confirm it using the OK key. Now the cursor will be flashing on the argon field ARGON, to confirm this operation activate the left arrow and the related field will be filled, and confirm it using the OK key, ARGON=SIM ATIVADO, repeat the previous step to program the quick cooling field RESFR=NÃO and also the auto start field AUTOS=NÃO. When these fields are filled, the selected function will be enabled and it will make part of the recipe. EX-9.

STATUS		REC°N	TEMP °C	VÁCUO	
T	E	L	A	4	0 0 1
P	R	É	Q	=	0 5 : 0 0
P	R	E	N	=	0 0 : 0 0
N	O	M	E	:	P A D R Á O . C E R Â M
ELEVADOR		MENSAGENS			

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7-11) Activate the PROGRAMA key again and a new screen TELA-5 will be visualized, in which the cursor will be flashing on the “D” field of the ascent function SOBE, move the cursor to the position (5) using the left and right arrows and confirm the operation by using the OK key. This option confirmation will be visualized in parenthesis on the ascent field SOBE. Repeat the previous step to program the descent function of the elevator DESCE, move the cursor to the position (3) using the left and right arrows and confirm the operation by using the OK key. EX-10.

STATUS		REC°N	TEMP °C	VÁCUO	
T	E	L	A	5	0 0 1
S	O	B	E	(5)	E L E V A D O R
D	E	S	C	E (3)	M D 1 2 3 4 5 6 7 8 9
N	O	M	E	:	M D 1 2 3 4 5 6 7 8 9
ELEVADOR		MENSAGENS			

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
7-12) Activate the PROGRAMA key again and the system will be back to the programming beginning TELA-1. At this stage the recipe is ready to be saved GRAVADA in the memory, check the step-by-step for saving recipes PG-24-item-11. To send the recipe to the platform just activate the monitoring key MONITORAÇÃO.

8 Programming pressed ceramic recipes

Ex-2 recipe model – pressed ceramic

EXEMPLO DE RECEITA PARA CERÂMICA CONVENCIONAL										
TEMP-1 °C	VEL-1 Fixo	PAT-1 Elevador	TEMP-2 °C	VEL-2 °C/min	PAT-2 min	VAC-2	TEMP-3 °C	VEL-3 °C/min	PAT-3 min	VAC-3 min
800	60	BASE	960	55	29:00	29:00	0	0	00:00	0

PRÉ. QUEIMA min	PRENSAGEM min	ARGON		RESFRIAMENTO		AUTO-START		ELEVATOR	
		Campo preenchido	Campo preenchido	Campo preenchido	Campo preenchido	UP	DOWN		
00:00	09:00		NO	X	YES		NO	D	D

8-1) Start the programming by reading the standard ceramic recipes (01); it will serve as reference for programming the other recipes.

8-2) Activate the LER (read) key and choose the type of work to be performed, ceramic or sintering. To alternate between the two options use the OK key and confirm the operation again using the LER (read) key. Ex4.

STATUS	RECºN	TEMP °C	VÁCUO	ELEVADOR	MENSAGENS	LEITURA DAS RECEITAS DE CERÂMICA	EX - 4
L E R							
0 1 à 7 9	0 1 à 7 9	C E R Â M I C A					

8-3) On the next screen choose receita 01 typing número 1 on the REC field, or navigate through the options list using the esquerda (left) or direita (right) arrows, the name of the recipe will be shown on the messages field (PADRÃO CERAM) EX-5, at this time you should confirm the option by activating the LER (read) key again to complete with the operation. The recipe will automatically be sent to the platform and when there is no recipe recorded, the message (RECEITA VAZIA) will be on the display.

10-3) On the next screen choose recipe 02 typing number 02 on the REC field or navigate through the options list using the left or right arrows, the recipe name will be on the message field (. . . OPACO . . .). At this stage, confirm the operation activating the LER key again to compete it. The recipe will be automatically sent to the furnace platform. EX-27.

STATUS	RECºN	TEMP °C	VÁCUO	LEITURA CERÂMICA	EX - 27
L E R	R E C : 0 0 2	C E R Â M I C A			
N O M E :	. . . O P A C O . . .				
ELEVADOR	MENSAGENS				

10-4) After reading the recipe the system will return to the monitoring screen. EX-30

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10-5) EX-28- Pressed ceramic burn. Activate the LER (read) key and choose the type of work to be performed, ceramic (cerâmica) or sintering (sinterização). To alternate between the two options use the OK key and confirm the operation again using the LER (read) key. Ex4

STATUS	RECºN	TEMP °C	VÁCUO	LEITURA DAS RECEITAS DE CERÂMICA	EX - 4
L E R					
0 1 à 7 9	0 1 à 7 9	C E R Â M I C A			
8 0 à 9 9	8 0 à 9 9	S I N T E R I Z A Ç Ã O			
ELEVADOR	MENSAGENS				

10-6) On the next screen choose recipe 03 typing number 03 on the REC field or navigate through the options list using the left or right arrows, the recipe name will be on the message field (. . . PRENSADA . . .) At this stage, confirm the operation activating the LER key again to compete it. The recipe will be automatically sent to the furnace platform. EX-29.

Pressed ceramic burn. REC-03

STATUS	RECºN	TEMP °C	VÁCUO	LEITURA CERÂMICA	EX - 29
L E R	R E C : 0 0 3	C E R Â M I C A			
N O M E :	. . . P R E N S A D A . . .				
ELEVADOR	MENSAGENS				

10-7) After reading the recipe the system will return to the monitoring screen. EX-31

STATUS		REC°N	TEMP °C		VÁCUO	
T	E	L	A	3	0 8 0	
		T	3	= 0 0 0 0	P A T 3 = H 0 : 0 0	
V	E	L	3	= 0 0 . 0		
N	O	M	E :		P A D R Ã O . S I N T E R .	
ELEVADOR		MENSAGENS				

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9-9) Activate the PROGRAMA key again and a new screen TELA-4 will be visualized, in which the cursor will be flashing on the "D" field of the ascent function SOBE, if necessary, move the cursor to the desired position using the left and right arrows and confirm the operation by using the OK key. This option confirmation will be visualized in parenthesis on the ascent field SOBE. Repeat the previous step to program the descent function of the elevator DESCE, move the cursor to the position (2) using the left and right arrows and confirm the operation by using the OK key. EX-24.

STATUS		REC°N	TEMP °C		VÁCUO	
T	E	L	A	4	0 8 0	
S	O	B	E	(D)	M D 1 2 3 4 5 6 7 8 9	
D	E	S	C	E (2)	M D 1 2 3 4 5 6 7 8 9	
N	O	M	E :		P A D R Ã O . S I N T E R .	
ELEVADOR		MENSAGENS				

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9-10) Activate the PROGRAMA key again and the system will be back to the programming beginning TELA-1. At this stage the recipe is ready to be saved GRAVADA in the memory, check the step-by-step for saving recipes PG-24-item-11. To send the recipe to the platform just activate the monitoring key MONITORAÇÃO.

10 Reading recipes

10-1) Reading recipes from the database.

Ex-27 – Conventional ceramic opaque burn. REC-02

Ex-28 – Pressed ceramic burn. REC-03

Ex-32 – Sintering and infiltration burn. REC-81

10-2) EX-27 - Opaque burn. Activate the LER (read) key and choose the type of work to be performed, ceramic (cerâmica) or sintering (sinterização). To alternate between the two options use the OK key and confirm the operation again using the LER (read) key. Ex4

STATUS		REC°N	TEMP °C		VÁCUO	
L	E	R				
0	1	à	7	9		
8	0	à	9	9		
					C E R Â M I C A	
					S I N T E R I Z A Ç Ã O	
ELEVADOR		MENSAGENS				

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8-4) After the standard recipe is read the system will return to its monitoring screen. EX- 11.

STATUS		REC°N	TEMP °C		VÁCUO	
M	O	N	I	T	0 0 1	0 6 0 0
T	1	=	0	6 0 0		
V	E	L	1	= 0 6 0		
S	D	.	D	D	P A D R Ã O . C E R Â M	
ELEVADOR		MENSAGENS				

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8-5) To start programming a new recipe you should activate the PROGRAMA key, and the display will show the message TELA1, the cursor will be flashing over the field related to temperature T1=600. At this time we have two options, we can simply type new values on the already filled fields or we erase all values using the LIMPA REC key, then we can type on the empty fields. EX-12.

STATUS		REC°N	TEMP °C		VÁCUO	
T	E	L	A	1	0 0 1	
T	1	=	0	0 0 0		
N	O	M	E :		
ELEVADOR		MENSAGENS				

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Note: For PRESSED CERAMIC EX-2, we will use the empty fields activating the function LIMPA REC.

8-6) On the screen TELA-1, now with the empty fields, you will see the cursor flashing on the T1=0000 field. Type the new value T1=0800 following the steps presented for standard pressed ceramic recipe EX-13.

STATUS		REC°N	TEMP °C		VÁCUO	
T	E	L	A	1	0 0 1	
T	1	=	0	8 0 0		
N	O	M	E :		
ELEVADOR		MENSAGENS				

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8-7) Activate the PROGRAMA key again and a new screen TELA-2 can be visualized, in which the cursor will be flashing on the T2=000 field. Type the new

value of T2=960 and confirm using the OK key, then the cursor will be flashing on the velocity field VEL2=000, type the new value of VEL2=55 and confirm it using the OK key. Now type the new value for the baseline time PAT2=29:00 and vacuum time VAC2=29:00, these times are programmed using minutes. Confirm the operation by using the OK key. EX-14.

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8-9) Activate the PROGRAMA key again and a new screen TELA-4 will be visualized, in which the cursor will be flashing on the pre burn field PRÉQ=00:00, keep it at ZERO and confirm it using the OK key, then the cursor will be flashing on the press field PREN=00:00, type the new value of PREN=09:00 and confirm it using the OK key. Now the cursor will be flashing on the argon field ARGON, keep it deactivated ARGON=NÃO and confirm this operation using the OK key. To program the quick cooling field RESFR, activate the left arrow to fill in the field RESF=SIM and confirm this operation using the OK key. The cursor will then be flashing on the auto start field, keep it deactivated AUTOS=NÃO and confirm it by using the OK key. When these fields are filled, the selected function will be enabled and it will make part of the recipe. EX-16.

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8-10) Activate the PROGRAMA key again and a new screen TELA-5 will be visualized, in which the cursor will be flashing on the "D" field of the ascent function

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9-6) On the screen TELA-1, now with the empty fields, you will see the cursor flashing on the T1=0000 field. Type the new value T1=220 and confirm using the OK key, then the cursor will be flashing on the velocity field VE1=00.0, type the new value of VEL1=07.0 and confirm it using the OK key. Now type the new value for the baseline time PAT1=H0:10 this time is programmed using hours. Confirm the operation by using the OK key. EX-21. Follow the step-by-step of the standard sintering recipe EX-3

Note: The sintering velocities can be fractional, values lower than 1°C/minute, minimum 0.1°C/minute.

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9-7) Activate the PROGRAMA key again and a new screen TELA-2 can be visualized, in which the cursor will be flashing on the T2=000 field. Type the new value of T2=1100 and confirm using the OK key, then the cursor will be flashing on the velocity field VEL2=00.0, type the new value of VEL2=40.0 and confirm it using the OK key. Now type the new value for the baseline time PAT2=H1:00, this time is programmed using hours with minutes visualization. Confirm the operation by using the OK key. EX-22.

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9-8) Activate the PROGRAMA key again and a new screen TELA-3 can be visualized, in which the cursor will be flashing on the T3=000 field. For this recipe, stage three (três) will not be used, so all the fields on this screen should be at ZERO EX – 15.

T3=000, VEL3=000, PAT3=00:00, VAC3=00:00
Confirm the operation by using the OK key. EX-23.

9-2) Activate the LER (read) key and choose the type of work to be performed, ceramic or sintering. To alternate between the two options use the OK key and confirm the operation again using the LER (read) key. Ex-19

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EX - 19

9-3) On the next screen choose receita 80 typing número 80 on the REC field, or navigate through the options list using the esquerda (left) or direita (right) arrows, the name of the recipe will be shown on the messages field (PADRÃO SINTER) EX-20, at this time you should confirm the option by activating the LER (read) key again to complete with the operation. The recipe will automatically be sent to the platform and when there is no recipe recorded, the message (RECEITA VAZIA) will be on the display. It is not possible to read (LER).

STATUS		REC°N	TEMP °C	VÁCUO		
L	E	R				
R	E	C :	0	8	0	S I N T E R I Z A Ç Ã O
A	L	F	A			
N	O	M	E :	P	A	D
ELEVADOR		MENSAGENS				

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EX - 20

9-4) After the standard sintering recipe is read the system will return to its monitoring screen. EX- 25

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9-5) To start programming a new recipe you should activate the PROGRAMA key, and the display will show the message TELA1, the cursor will be flashing over the field related to temperature T1=600. At this time we have two options, we can simply type new values on the already filled fields or we erase all values using the LIMPA REC key, then we can type on the empty fields. EX-26.

Note: For SINTERING EX-3, we will use the empty fields activating the function LIMPA REC.

SOBE, if necessary, move the cursor to the desired position using the left and right arrows and confirm the operation by using the OK key. This option confirmation will be visualized in parenthesis on the ascent field SOBE. Repeat the previous step to program the descent function of the elevator DESCE, move the cursor to the position (D) using the left and right arrows and confirm the operation by using the OK key. EX-17.

STATUS		REC°N	TEMP °C	VÁCUO			PROGRAMAÇÃO
T	E	L	A	5	0	0	1
S	O	B	E	(D)	M	D	1
D	E	S	C	E	(D)	M	D
N	O	M	E	: . . .	P	R	E
ELEVADOR		MENSAGENS					SEQUENCIAL DE CERÂMICA PRENSADA


EX - 17

8-11)Activate the PROGRAMA key again and the system will be back to the programming beginning TELA-1. At this stage the recipe is ready to be saved GRAVADA in the memory, check the step-by-step for saving recipes PG-24-item-11. To send the recipe to the platform just activate the monitoring key MONITORAÇÃO

9 Programming sintering recipes.

Ex-3 recope model – sintering/infiltration

EXEMPLO DE RECEITA PARA SINTERIZAÇÃO									
TEMP-1 °C	VEL-1 °C/min	PAT-1 H/min	TEMP-2 °C	VEL-2 °C/min	PAT-2 H/min	TEMP-3 °C	VEL-3 °C/min	PAT-3 H/min	ELEVATOR
200	7.0	00:10	1100	40.0	01:00	0	0	00:00	D 2

9-1) Start the programming reading the standard sintering recipe Receita (80), it will serve as reference for programming the other recipes.